

Infrared forehead thermometer

Item No. 00164


PZN (D) 16833831

Description

Contactless infrared forehead thermometer with fever alarm and memory function - fast, accurate and easy to use. Ideal for measuring the body temperature of adults, children and infants. It can be used indoors and outdoors and is also suitable for measuring object temperatures.

The thermometer is equipped with a backlight for use in the dark and a convenient auto shut-off feature. The color-coded display with its three color warning levels makes it easy to interpret temperature readings.

- ✓ Measurement of body and object temperatures
- ✓ Fast measurement in just one second
- ✓ Three color warning levels (green, orange, red)
- ✓ Automatic switch-off after 60 seconds



Technical data

Display	LCD digital
Measuring points	Forehead and object surface
Measurement ranges	Body Mode: 32,0°C – 43,0°C (89,6°F – 109,4°F) Object mode: 0,0°C – 100,0°C (32,0°F – 212,0°F)
Temperature unit	°C / °F
Measurement unit	0,1°C / 0,1°F
Measurement duration	1 Second
Memory function	Last measured value
Power source	2 pcs AAA alkaline batteries
Automatic power off	After 60 seconds
Weight	ca. 95 g (batteries not included)
Dimensions	95 x 45 x 155 mm
Battery life	ca. 500 temperature measurements
Content	Thermometer, operating instructions, protective bag, battery

